




MAKE EVERY MINUTE COUNT

COMPLETE PATHOLOGY SOLUTIONS



all the
right tools.

always ready
for you.

Discover the next breakthrough — and see the impact you can make on the world.

With our complete pathology solution, you get the right tools, products and superior customer support throughout your entire workflow. Now you can have confidence in the quality of every result and look boldly toward tomorrow.

**With our one-stop pathology solution, there's
no limit to what you can do:**

find the diagnosis faster, save a life, change the world.



Sample Collection and Handling

- Dissection Tools
- Ancillary Components (slides, coverslips, etc.)



Immunohistochemistry

- SNAP i.d.® Protein Detection System
- Chromogenic Substrates
- Antigen Retrieval
- Research Use Antibodies
- Cell Marque™ IVD Antibodies



Fixation and Processing

- Fixatives
- Solvents
- Mounting Media
- Embedding Media



Water Purification Solutions



Tissue Staining

- Routine Stains
- Special Stains
- TISSUE-TROL™ Control Slides
- Automated Stainer

Single products used as Fixatives, Tissue processing agents, Mounting Media and Embedding Media and many Stains belong to the group of *In Vitro* Diagnostic (IVD) medical devices. According to US FDA regulation they are classified being Class I, exempt IVD's. Most of them comply with the EU IVD Directive, bearing the CE logo on the label.

Products included in this brochure are not available in all countries.
Contact your local sales representative or dealer for details.

Sample Collection and Handling

Dissection Tools

A selection of dissection instruments with consistent product quality you can trust.

- **Scissors** – Selection of stainless steel Micro-dissection, Mayo and Sharp forms all available with straight or curved blades
- **Scalpels** – Choose the system that suits you best, disposable with plastic handles or disposable stainless steel blades for use with reusable handles.
- **Forceps and Tweezers** – There is a wide variety of point styles and designs including Dumont, hemostatic and mosquito.



Cat. No.	Product Description
Z162507	All-purpose tweezer set, 5 tweezers in a vinyl pouch, pkg of 1 set
F2394	Bayonet forceps, stainless steel
Z692395	Disposable scalpel No 10, individually wrapped
Z675776	DissecTable™ reversible dissecting board polyethylene, pkg of 1 ea
Z265993	Dissecting scissors, Mayo, curved
Z265977	Dissecting scissors, Sharp, straight
F6521	Jewelers forceps, Dumont No. 5
F3767	Micro-dissecting forceps, very fine point
S3146	Micro-dissecting scissors, stainless steel straight, sharp point
S2646	Scalpel blades, No. 10
S2771	Scalpel blades, No. 11
S2896	Scalpel handles, No. 3, stainless steel
Z168866	Scissor-handle forceps, Kelly hemostatic, straight
Z509140	Scissor-handle forceps, straight, Occulating
Z225304	Stainless steel forceps, High-precision
Z168793	Stainless steel forceps, Dissecting
Z168785	Stainless steel forceps, Fine tip, curved
Z168777	Stainless steel forceps, Fine tip, straight
Z680214	Titanium tweezers, Style 2A
T5415	Tweezers, Style #00, sharp, anti-magnetic stainless steel
T5790	Tweezers, Style #4, needle-sharp, anti-magnetic stainless steel
T4662	Tweezers, Style #5, Fine needle-sharp, Inox alloy
T4537	Tweezers, Style #5, Fine needle-sharp, non-magnetic stainless steel
T4412	Tweezers, Style #5, Fine needle-sharp, non-magnetic alloy

Successfully start your pathology research with our one-stop product resources at your side. Now you have even more support to discover the next breakthrough. Find new ways to identify, control and prevent disease—and ultimately, improve patient care.

Collection Vessels

To make sample collection as simple as possible we have a range of plain vessels in a wide variety of volumes and designs. To supplement this we have tamper evident cap options and a selection of specialist vessels prefilled with preservatives.

Cat. No.	Product Description
Z740425	HistoTainer™ I prefilled specimen container, tamper evident, volume 20 mL
Z740427	HistoTainer™ I prefilled specimen container, tamper evident, volume 60 mL
Z186821	All-purpose sampling containers capacity 4 1/2 oz
Z186848	All-purpose sampling containers capacity 8 oz
Z645346	Bijou sample container plain label
BR61815	BRAND® jar with push-on lid polypropylene, capacity 40 mL
BR61820	BRAND® jar with push-on lid polypropylene, capacity 80 mL
BR61960	BRAND® jar with screw cap low-density polyethylene, capacity 30 mL
BR61950	BRAND® jar with screw cap low-density polyethylene, capacity 5 mL
BR61970	BRAND® jar with screw cap low-density polyethylene, capacity 90 mL
BR62310	BRAND® jar with screw cap polypropylene, capacity 30 mL
BR62315	BRAND® jar with snap-on PE-LD lid polystyrene, capacity 25 mL
BR62020	BRAND® PFA sample jar with screw cap volume 180 mL
BR62005	BRAND® PFA sample jar with screw cap volume 30 mL
BR62015	BRAND® PFA sample jar with screw cap volume 90 mL
Z740379	Eco-Friendly SpecTainer™, tamper evident, biodegradable, size 120 mL, sterile
Z740376	Eco-Friendly SpecTainer™, tamper evident, biodegradable, size 60 mL, non-sterile
BAF175750003	Snap Cap Vials capacity 11 mL
BAF175750040	Snap Cap Vials capacity 148 mL
BAF175750013	Snap Cap Vials capacity 48 mL
Z740422	Tamper proof specimen containers, size 1000 mL
Z740420	Tamper proof specimen containers, size 300 mL
Z740424	Tamper proof specimen containers, size 5700 mL
Z645362	Universal sample container with flow seal cap plain label



Discover our complete product listing at [SigmaAldrich.com/histolabware](https://www.sigmaaldrich.com/histolabware)

Slides/Coverslips

Explore our comprehensive range of simple slides and coverslips along with cavity slides and specialist coatings.

- **Microscope Slides** – 4 sizes, 7 end options, 3 colors and available with cut or ground edges
- **Coverslips** – 9 sizes in 2 thicknesses
- **Specialist Options** – Cell chamber slides, coated slides, stick on chambers for microscope slides are just a few examples



Cat. No.	Product Description
BR475505	BRAND® cavity slides, 76 mm × 26 mm, with 1 concavity
BR470820	BRAND® cover glass, 24 mm × 60 mm
BR474702	BRAND® microscope slide, cut edges, twin frosted end
BR474743	BRAND® microscope slide, ground edges
Z681164	Colored Microscope slides, blue
CLS40042	Corning® epoxide coated slides
CLS40004	Corning® GAPST™ II coated slides
CLS294775X25	Corning® microscope slides, plain, 75 mm × 25 mm
CLS294775X50	Corning® microscope slides, plain, 75 mm × 50 mm
C9802	Cover glass, 22 mm × 22 mm
C9056	Cover glass, 24 mm × 60 mm
GBL622305	Grace Bio-Labs CoverWell™ perfusion chambers, 8 chambers 9 mm diam.
GBL716022	Grace Bio-Labs HybriSlip™ hybridization cover, 22 mm × 60 mm
GBL666504	Grace Bio-Labs Press-To-Seal silicone isolator, 20 mm diam.
GBL654004	Grace Bio-Labs SecureSeal™ imaging spacer, 1 well, 13 mm diam.
Z688533	Microscope slide-grids, disposable
S8400	Microscope slides, frosted one end, 25 mm × 75 mm
S8902	Microscope slides, plain, 25 mm × 75 mm
C7057	Nunc® Lab-Tek® II Chamber Slide™ system
Z734756	Nunc® Lab-Tek® II chambered coverglass, 4 wells, 1.7 cm ² /well
Z734853	Nunc® Lab-Tek® II chambered coverglass, 8 wells, 0.7 cm ² /well
Z692255	Superfrost microscope slides, white end, 76 mm × 26 mm


Staining Labware

We offer a variety of staining dishes and Coplin jars suitable for pathology, microscopy, microbiology and hematology. Dishes, covers and racks are manufactured from PMP, Black POM or soda-lime glass and resist stains from eosin or hematoxylin. Staining dishes are available in many configurations to accommodate up to 50 slides, with various plastic, glass and stainless steel racks.

For immunostaining, we have a selection of staining trays that maintain humidity and reduce evaporation as well as PAP pens to ensure the minimum volume of precious reagent is used on each slide.

We can even help you if you end up staining your hands rather than your samples with a very effective stain remover.

Cat. No.	Product Description
BR472800	BRAND® glass staining trough, Coplin pattern, soda-lime glass
BR472700	BRAND® glass staining trough, Hellendahl extended, soda-lime glass
BR472500	BRAND® glass staining trough, Schiefferdecker pattern, soda-lime glass
BR474400	BRAND® plastic staining trough, Hellendahl pattern, with lid
BR472000	BRAND® staining tray for glass staining trough, soda-lime glass
BR471400	BRAND® staining tray, black, size 76 mm x 26 mm, for 25 slides, POM
BR471500	BRAND® staining trough, black
BR471800	BRAND® staining trough, complete with tray, black, for 25 slides
BR474305	BRAND® transparent plastic staining tray, for 20 slides, use with transparent plastic staining trough
Z675776	DissecTable™ reversible dissecting board polyethylene, pkg of 1 ea
Z151564	Erado-Sol® liquid stain remover
Z358290	Nalgene® staining box, L x W x H 12.5 cm x 12.5 cm x 5 cm
Z358304	Nalgene® staining box, L x W x H 22.5 cm x 22.5 cm x 5 cm
Z377821	PAP pen for immunostaining, 5 mm tip width
S5641	Scienceware® Coplin staining jar with screw cap, high-density polyethylene
Z674664	Slide spinner for slide drying
Z670146	StainTray with black lid
Z670138	StainTray with clear lid
S6141	Wheaton Coplin staining jars, 20 slide staining dishes complete with covers, slide rack and handle
S5516	Wheaton Coplin staining jars, 5 slide jar with ground glass cover
S5891	Wheaton Coplin staining jars, 5 slide jar with plastic screw cap
S6016	Wheaton Coplin staining jars, 8 slide jar with glass cover
Z103969	Wheaton staining dishes includes: dish, rack, and handle
S5017	Wheaton staining dishes, 20 slides, rack only
S4642	Wheaton staining dishes, for 20 slides, dishes only
S4892	Wheaton staining dishes, for 20 slides, dishes with covers only



Fixation and Processing



Fixatives

Our fixatives are pre-formulated for direct use on tissue specimens allowing for efficient and safe protocols.

Cat. No.	Product Description
HT101128-4L HT10132-1L	Bouin's solution – Prepared with saturated picric acid, formaldehyde and acetic acid. Excellent fixative for preserving soft and delicate structures. Used as a mordant in various trichrome procedures.
252549-25ML 252549-500ML 252549-6X500ML 252549-1L 252549-4L 252549-4X100ML 252549-100ML	Formaldehyde solution, ACS reagent, 37 wt. % in H ₂ O – This product is a solution of approximately 37% by weight of formaldehyde gas in water. The product contains 10–15% methanol as a stabilizer to prevent polymerization.
HT5011-48X15ML HT5012-24X60ML HT5014-24X120ML HT501128-4L HT501320-9.5L HT501640-19L	Formalin solution, neutral buffered, 10% – Case containers are polypropylene with white leak-proof lids and are partially pre-filled with 10% NBF. Case sizes are 48×15 mL, 24×60 mL and 24×120 mL. Popular general purpose histological fixative.
A5472-1GA	Formalin Free Tissue Fixative 🌱 – A less toxic alternative to formalin. Performs well with PCR, <i>in situ</i> hybridization and immunohistological staining.
F5554-4L F5554-19L	Formalin, 10%, Neutral Buffered, Wintergreen Scented – Intended for use as a general purpose histological fixative. Wintergreen oil has been added to the standard formulation to provide a wintergreen scent.
F5304-4L	Formalin, 10%, Neutral Buffered with 0.03% Eosin – Intended for use as general histological fixative. The addition of the eosin to the formulation allows the user to more easily visualize small fragments of tissue during the transfer of the tissue from the fixative to the tissue cassette.
103999	Formaldehyde solution min. 37% free from acid stabilized with about 10% methanol and calcium carbonate for histology, 1L, 2.5L, 25L
65115	HARLECO® Formalin, Carsons. 1:10 dilution of neutral buffered formalin used for tissue fixation in histology. Can be used for tissue specimens in electron microscopy, 20L
H2904-100ML H2904-1L	HistoChoice® Tissue Fixative 🌱 – Designed for the molecular biologist. Specially formulated to preserve antigenic sites for antibody probes and nucleic acid sites for <i>in situ</i> hybridizations. Fixed sections do not require pre-digestion or other recovery procedures to make important sites available.
P6148-500G P6148-1KG P6148-5KG	Paraformaldehyde, reagent grade, crystalline – Recommended for use in the preparation of formalin fixatives for tissues or cells when the samples are to be used in fluorescence studies.
101736	OSTEOMOLL® fixation and decalcification in one step for bone and hard tissue: 1L, 2.5L
100974	Ethanol absolute, for analysis EMSURE® ACS, ISO, Reag. Ph Eur: 1L, 5L
P6744-1GA	Picric acid solution, 1.3% in H ₂ O (saturated) – For use as a colorimetric reagent in creatinine determination assay. For use as a general cytoplasmic stain, connective tissue fiber stain, and as a fixative for tissues.
Z2902-3.75L	Zinc Formalin Fixative – Excellent morphological preservation of nuclear and cytoplasmic components. Can replace neutral buffered formalin for routine tissue and immunohistochemical procedures. It is a non-precipitating fixative which can be used with automated and manual methods.

New Line of Reagent Alcohol (RGA)

MilliporeSigma is proud to introduce a new line of ethanol products. We now offer Tax-free Pure Ethanol and Reagent Alcohol (RGA) in a wide range of packaging options for your histology and pathology applications.

This new line of products is supported by a robust distribution network to ensure expedient delivery and low cost of delivery.

- Range of sizes and prices to suit any budget
- Dense distribution network – low shipping cost
- Not subject to federal excise tax
- Requires letterhead approval from the TTB for volumes greater than 4 liters

Reagent Alcohol (RGA)

Cat. No.	Product Description
793175	100% (200 proof), pkg sizes: 4x1GA, 5GA, 55GA, 270 GA
793183	95% (190 proof), pkg sizes: 4x1GA, 5GA, 55GA, 270 GA
793191	80% (160 proof), pkg sizes: 4x1GA
793213	70% (140 proof), pkg sizes: 4x1GA
752X	HARLECO® Acid Alcohol, 70% (3% HCl), 1L
753X	HARLECO® Acid Alcohol, 95% (5% HCl), 1L

Pure Ethanol

Cat. No.	Product Description
792799	190 proof, ACS reagent, meets USP testing specifications, pkg sizes: 24x1 PT, 20x1 GA, 120x1 PT, 4x1 GA, 55 GA, 5 GA
792780	200 proof, ACS reagent, meets USP testing specifications, pkg sizes: 24x1PT, 4x1GA, 5GA, 55GA, 120 x 1Pt , 20 x 1 GA

Alternative Tissue Clearing Agents

HistoChoice® Clearing Agent

HistoChoice® is a patented, complex mixture of agents that allow tissues to retain their structure, antigenic and nucleic acid binding sites. This formulation is non-citrus, and contains no mercury, formaldehyde or glutaraldehyde.

HistoChoice® Clearing Agent is suitable as an alternative to toluene and xylene for dewaxing paraffin tissue sections. Prepared slides are suitable for antibody probing applications and in situ hybridizations.

Features

- Non-toxic and odorless
- Preserves native antigenic and nucleic acid sites
- Tissues retain original structure
- Replaces common fixatives that are often toxic

Cat. No.	Product Description
H2779-1L	HistoChoice® clearing agent



Xylene Substitute

Non-toxic xylene replacement. Compatible with all tissue processors. May not be compatible with all permanent mounting medias.

Cat. No.	Product Description
A5597-1GAL	Xylene substitute
A5597-4X1GAL	
109843	Neo-Clear® is an aliphatic hydrocarbon clearing agent. A green substitute for xylene: 5L, 4 x 5L, 25L

Classic Tissue Clearing Agents

Cat. No.	Product Description
534056	Xylenes, histological grade
472476	Chloroform, ACS reagent, ≥99.8%(GC)
179418	Toluene, ACS reagent, ≥99.5%(GC)
M29908	Methyl benzoate, ≥ 99% (GC)
M6752	Methyl salicylate, ≥99% (GC)

Discover our complete product listing at SigmaAldrich.com/ethanol


Tissue Processing

Cat. No.	Product Description
101736	OSTEOMOLL® fixation and decalcification in one step for bone and hard tissue: 1L, 2.5L
101728	OSTEOSOFT® Decalcifying solution for sensitive hard material and bone marrow biopsies. Can be used for specimens in immunohistochemistry: 1L, 10L
65039	Mollifex™ Tissue Softener. Softening agent for use with paraffin embedded tissues. Ideal for hard to section blocks and brittle tissues: 500mL

Mounting and Embedding Labware

Our reagent offering is supported by a wide range of assorted labware to aid in the embedding and mounting your samples.

Cat. No.	Product Description
BR112805	BRAND® forcep with pointed ends PMP, L 115 mm
BR113105	BRAND® forcep with pointed ends polyoxymethylene, L 100 mm
BR113125	BRAND® forcep with rectangular ends PTFE, L 100 mm
E4140	Embedding molds 19 Flat bullet-shaped unindexed cavities, clear
E4390	Embedding molds 19 Flat bullet-shaped unindexed cavities, maroon
E4015	Embedding molds 26 Flat bullet-shaped cavities double-indexed with numbers and letters, clear
E4265	Embedding molds 26 Flat bullet-shaped cavities double-indexed with numbers and letters, maroon
G2528	Glass strips for thick section (Histology) knives
G2403	Glass strips for thin section (EM) knives
Z648183	Grid coating pen provides a thin layer of adhesive
Z676225	Micro-cutting tool set for microscopy scale applications
Z374415	Scienceware® microscope slide dispenser Holds approximately 72 slides
Z374407	Scienceware® microscope slide holder accepts up to 10 single or double slides standing on end
Z374385	Scienceware® slide box, hinged provides the utmost protection for up to 100 - 1×3 in. slides
Z648205	Securline lab marker II lab marker, for use on plastic embedding cassettes and frosted end of slides
Z674664	Slide spinner for slide drying AC input 110 V
Z648191	STATMARK permanent marker for cassette & slide marking
Z648175	Tissue capture pen one pen for two applications
Z672122	Tissue processing/embedding cassettes with lid



Discover our complete product listing at SigmaAldrich.com/histolabware

Mounting Media

We have a range of natural and synthetic mounting media that provide an ideal platform to mount your valuable samples.


Cat. No.	Product Description
I1161-30ML	ImmunoHistoMount™ - A permanent aqueous mounting medium designed for tissue sections and cell smears with peroxidase and alkaline phosphatase chromogens that cannot be dehydrated with organic solvents.
C1795-25ML C1795-100ML	Canada Balsam - A natural mounting medium obtained from balsam fir trees. The optical properties are nearly identical with those of glass.
C9368-30ML	CC/Mount™ - An aqueous based permanent mounting medium, similar to Crystal Mount™, for the permanent preservation of tissue sections stained with peroxidase and alkaline phosphatase based systems as well as with various fluorescent dyes.
GG1-15ML GG1-10X15ML	Glycerol Gelatin - Aqueous slide mounting media for histological use. Product must be warmed to 55-60 °C before use; at ~37 °C the product returns to semi-solid state.
O8015-30ML O8015-100ML	Organo/Limonene Mount™ 🌿 - This ready-to-use mounting medium is made with limonene, a natural product from orange peels. It is good for coverslipping stained tissues and cell smears that can be dehydrated with organic solvents. This mounting medium is suitable for organic solvent resistant immunohistochemistry (IHC) chromogens as well as an excellent choice for mounting H and E stained slides.
64969	HARLECO® Krystalon™ Mounting Medium: 300mL, 500mL
109016	Neo-Mount™ is a non-aqueous, user-friendly mounting medium designed to work with Neo-Clear® xylene substitute: 100mL, 500mL



GO GREENER HISTOLOGY & HEMATOLOGY

MilliporeSigma is committed to bringing safer greener alternatives to our customers.





Embedding Media

Cat. No.	Product Description
EM0100	JB-4 Embedding Kit – JB-4 is a water-soluble, GMA based, plastic resin kit intended for use in the preparation of embedded samples for high resolution light microscopy. Yields 0.5µ-2µ sections
EM0200	Osteo-Bed Bone Embedding Kit – Suitable for use with large and small mineralized (undecalcified) bone sections. Yields clear, hard blocks for cutting sections. Not water- soluble.
O8639	Osteo-Bed Bone Embedding Solvent – For use with Osteo-Bed Bone Embedding Kit. Solvent for removal of plastic from sections prior to rehydration and staining.
P0091	PolyFreeze – Water-soluble support matrix and a form of embedding medium for frozen sectioning.
EM0300	Spurr Low Viscosity Embedding Kit – Low viscosity for exceptional penetration of tissue and other difficult matrices.
SHH0023-6X120ML	PolyFreeze Tissue Freezing Medium – A fast freezing support matrix for frozen tissue sectioning, Green
SHH0024-6X120ML	PolyFreeze Tissue Freezing Medium – A fast freezing support matrix for frozen tissue sectioning, Yellow
SHH0025-6X120ML	PolyFreeze Tissue Freezing Medium – A fast freezing support matrix for frozen tissue sectioning, Blue
SHH0026-6X120ML	PolyFreeze Tissue Freezing Medium – A fast freezing support matrix for frozen tissue sectioning, Clear
P3558-1KG	Paraplast® embedding media – Produces minimal tissue compression and wrinkle-free sections, and cuts to 4 µm thickness with excellent ribbon continuity.
P3683-1KG	Paraplast® Plus embedding media – Regular Paraplast with approx. 0.8% DMSO added for enhanced ease of infiltration and sectioning for tissues difficult to process with regular Paraplast.
P3808-1KG	Paraplast® Xtra embedding media – A unique blend of highly purified paraffin and low M.W. polymers for lower temperature infiltration and better compression resistance.
111609	Histosec® pastilles. Paraffin enriched with polymers: 1kg, 4 x 2.5kg, 25kg
115161	Histosec® pastilles without DMSO. Paraffin enriched with polymers: 4 x 2.5kg, 25 kg

Engineered Histology Bath Fluid

3M™ Novec™ Non-Flammable Bath Fluids

Environmentally Friendly and Non-Flammable Bath Fluid

Traditional fluids used for freezing histology samples (including isopentane) are highly flammable and toxic. This has implications for worker safety, operating procedures, building code compliance and liability. Today's laboratories demand a safer alternative without compromising performance. We are pleased to offer 3M™ Novec™ engineered fluids in convenient package sizes to meet your laboratory needs.

- Non-flammable, low toxicity
- High density = reduced frost
- Inert to embedding medium
- Inert to tissue sample
- Less volatile = lower evaporation
- Low global warming potential (GWP)
- Low odor

Cat. No.	Product Description
SHH0001	Description: 3M™ Novec™ 7000 Engineered Fluid, 1-methoxyheptafluoropropane, is a non-flammable, low global warming potential (GWP) heat transfer fluid capable of reaching -120°C. It is also useful as a direct expansion refrigerant
SHH0002	3M™ Novec™ 7100 Engineered Fluid, methoxy-nonafluorobutane (C4F9OCH3), is a clear, colorless and low-odor fluid intended to replace ozone-depleting substances (ODSs) and compounds with high global warming potential (GWP) in many applications

Routine

Tissue Staining

Differential/Counter Stains

The combination of the nuclear stain hematoxylin and the plasma counterstain eosin Y is the most useful stain in the field of histotechnology. Various formulations of hematoxylin and eosin, prepared with dyes certified by the Biological Stain Commission, are supplied as ready-to-use solutions. Dyes, chemicals and general-purpose laboratory reagents are also available for customers wishing to prepare their own staining solutions. The histology product line includes reagents for *in vitro* diagnostic use for clinical laboratories.

Eosins

General purpose cytoplasmic stain highlighting in pink to red. Eosin Y stains protein in the cytoplasm due to its cell membrane permeability and attraction of the negatively charged dye to positively charged amino acid residues of proteins in the cytoplasm and in collagen and muscle fibers. Available in aqueous and alcoholic formulations.

Cat. No.	Product Description
HT110116-500ML HT110132-1L HT110180-2.5L HT1101128-4L	Eosin Y solution, alcoholic – Certified Eosin Y, 0.5% (w/v) in acidified 90% ethanol
HT110316-500ML HT110332-1L HT110380-2.5L HT1103128-4L	Eosin Y solution, alcoholic, with phloxine – Certified Eosin Y, 0.1% (w/v) and certified phloxine B, 0.1% (w/v) in acidified ethanol
HT110216-500ML HT110232-1L HT110280-2.5L HT1102128-4L	Eosin Y solution, aqueous – Certified Eosin Y, 0.5% (w/v) in water, not acidified
586X	HARLECO® Eosin Y Solution, Alcoholic Saturated: 1L
588X	HARLECO® Eosin Y Solution, Aqueous, 1% PDC: 1L, 10L

Hematoxylin

General purpose nuclear stain highlighting in dark blue to violet. Hematoxylin stains oligonucleotides in the nucleus due to its attraction and binding to the negatively charged oligonucleotide phosphate backbone.

Cat. No.	Product Description	Cat. No.	Product Description
GHS116-500ML GHS132-1L GHS1128-4L	Hematoxylin Solution, Gill No. 1 – Used as a progressive cytology stain, used with hematoxylin and eosin staining, 2 g/L certified hematoxylin	MHS1-100ML MHS16-500ML MHS32-1L MHS80-2.5L MHS128-4L	Hematoxylin Solution, Mayer's – General purpose, progressive type, used as a counterstain for procedures such as immunohistochemistry or laser microdissection, 1 g/L certified hematoxylin
GHS216-500ML GHS232-1L GHS280-2.5L GHS2128-4L	Hematoxylin Solution, Gill No. 2 – May be used for cytology and/or histology as a progressive or regressive stain depending on length of staining time, 4 g/L certified hematoxylin	109249	Mayer's Hemalum solution: 500mL, 1L, 2.5L
GHS316-500ML GHS332-1L GHS380-2.5L GHS3128-4L	Hematoxylin Solution, Gill No. 3 – May be used for cytology and/or histology as a progressive or regressive stain depending on length of staining time, 6 g/L certified hematoxylin	HT1079-1SET	Weigert's iron hematoxylin solution, set – Recommended in lengthy or acidic stain procedures, resists discoloration, 1% hematoxylin in ethanol and Ferric chloride in dilute HCl
105174	Hematoxylin solution modified acc. to Gill III; 500mL, 1L, 2.5L	S5134-6X100ML	Scott's Tap Water Substitute Concentrate, 10x – Used as a "blueing reagent" in hematoxylin and eosin staining procedures, MgSO ₄ • 7H ₂ O, 200 g/L + NaHCO ₃ , 20 g/L
HHS16-500ML HHS32-1L HHS80-2.5L HHS128-4L	Hematoxylin Solution, Harris Modified – General purpose, regressive type, 7 g/L certified hematoxylin	A3179-1L A3429-4L	Differentiation Solution – Acidified alcohol solution for the differentiation of regressive hematoxylin stains, Acid Alcohol

Special Stains


Special Staining Kits

Special stains are used to identify suspected pathogens or demonstrate specific cellular components that aid pathologists in the evaluation of disease states. We offer the most commonly utilized special stain products, with individual reagents available as prepared solutions or in convenient kit form. For customers wishing to prepare their own staining solutions we offer stains and dyes, as well as other chemicals and general purpose laboratory reagents required.

Amyloid Stain, Congo Red Kit

Cat. No.: HT60-1KT

Amyloid protein is detected in tissue with Congo red, a metachromatic stain. Staining is intensified by pretreatment of tissue with alkaline sodium chloride.



Components:

- Congo red solution (kit only), 500 mL
- Sodium chloride solution, alcoholic (kit only), 500 mL
- Sodium hydroxide solution (kit only), 2 × 12 mL

Silver Stain (Modified GMS) Kit

Cat. No.: HT100A-1KT

Silver stain is an aid in distinguishing fungi, basement membrane and opportunistic organisms such as *Pneumocystis carinii* in tissue sections or smears.



Components:

- Borax solution, 100 mL
- Gold chloride solution, 100 mL
- Periodic acid solution, 100 mL
- Silver methenamine reagent, 12 vials
- Sodium thiosulfate solution, 100 mL

Periodic Acid-Schiff Kit (PAS)

Cat. No.: 395B-1KT

Staining patterns of lymphocytes are helpful in making therapeutic decisions in established cases of lymphocytic leukemia. Cellular elements which may be demonstrated with the PAS procedure include glycogen, fungal walls, basement membrane, certain epithelial sulfomucins and sialomucins, neutral mucosubstances, colloid material of the thyroid and the pars intermedia of the pituitary.



Components:

- Hematoxylin solution, Gill No. 3 (kit only), 100 mL
- Periodic acid solution (395-1), 2 × 100 mL
- Schiff's reagent (kit only), 4 × 50 mL

Trichrome Stain (Masson) Kit

Cat. No.: HT15-1KT

Collagen and muscle fiber stain differentially when treated sequentially with Biebrich scarlet-acid fuchsin, PTA/PMA and aniline blue. Cytoplasm and muscle fibers stain red whereas collagen displays blue coloration.



Components:

- Aniline blue solution (kit only), 250 mL
- Biebrich scarlet-acid fuchsin solution (HT15-1), 250 mL
- Phosphomolybdic acid solution (HT15-3), 250 mL
- PPhosphotungstic acid solution (HT15-2), 250 mL

Iron Stain Kit

Cat. No.: HT20-1KT

Deposits of loosely bound iron in erythrocytes, bone marrow, or tissue stain an intense blue when treated with acid ferrocyanide; according to the well-known Prussian blue reaction.



Components:

- Hydrochloric acid solution, 250 mL
- Potassium ferrocyanide solution, 250 mL
- Pararosaniline solution, 30 mL

Reticulum Stain Kit

Cat. No.: HT102A-1KT

Intended to demonstrate reticular fibers in paraffin sections. Useful in the differential diagnosis of certain type of tumors such as carcinomas, sarcomas, and lymphosarcomas.




Components:

- Ferric ammonium sulfate solution, 500 mL
- Gold chloride solution, 100 mL
- Oxalic acid solution, 100 mL
- Potassium permanganate solution, 100 mL
- Silver nitrate solution, 50 mL
- Sodium hydroxide solution, 100 mL
- Sodium thiosulfate solution, 500 mL

DNA staining acc. to Feulgen Kit for quantitative determination of DNA

Cat. No.: 107907

The standardized Feulgen reaction is the method of choice for determining DNA in histological and cytological specimens. The significance of DNA content in diagnosis and prognosis in undifferentiated tumors is highly recognized.



Components:

- 5 M Hydrochloric acid, 2 x 250 mL
- Schiff's reagent, 2 x 250 mL
- Sodium metabisulfite concentrate, 250 mL

Elastica van Gieson's Kit for detection of elastic fibers

Cat. No.: 115974

Brilliant results are achieved for the visualization of fibers using resorcin-fuchsin solution according to Weigert's in association with van Gieson's picrofuchsin solution and nuclear staining with Weigert's iron hematoxylin solution.




Components:

- Weigert's solution A – Alcoholic hematoxylin solution, 500 mL
- Weigert's solution B – Hydrochloric acid iron(III)nitrate solution, 500 mL
- Elastin acc. to Weigert's resorcin fuchsin solution, 500 mL
- Picrofuchsin solution acc. to van Gieson

Masson-Goldner Kit for visualizing connective tissue

Cat. No.: 100485

This "Masson-Goldner staining kit - for the visualization of connective tissue with trichromic staining" is used for human-medical cell diagnosis and serves the purpose of the histological investigation of sample material of human origin.



Components:

- Azophloxine solution, 500 mL
- Tungstophosphoric acid orange G solution, 500 mL
- Light green SF solution, 500 mL
- Acetic acid 10%, 500 mL

Weigert's iron hematoxylin kit

Cat. No.: 115973

Following nuclear staining with Weigert's iron hematoxylin, components such as muscle, cytoplasm and erythrocytes are stained with azophloxin and Orange G solution. Connective tissue is then counterstained with Light green SF solution.



Components:

- Azophloxine solution, 500 mL
- Tungstophosphoric acid Orange G solution, 500 mL
- Light green SF solution, 500 mL
- Acetic acid 10%, 500 mL

NEW

Fontana-Masson Stain Kit



Cat. No.: HT200-1KT

Fontana-Masson staining is used for the demonstration of argentaffin substances such as melanin, argentaffin granules of carcinoid tumors and some neurosecretory granules. The 30X faster microwave procedure reduces incubation time: melanin stain time is reduced to 1 minute, argentaffin stain time is reduced to 2 minutes.



Components:

- Silver Nitrate Solution, 50 mL
- Gold Chloride Solution, 125 mL
- Sodium Thiosulfate Solution, 125 mL

30x Faster!

Silver Plating Kit acc. to von Kossa Kit for the detection of calcium deposits

Cat. No.: 100362

This "Silver plating kit acc. to von Kossa - for detection of microcalcification" is used for human-medical cell diagnosis and serves the purpose of the histological investigation of sample material of human origin.



Components:

- Staining kit contains, 2 x 100-ml bottles
- Silver nitrate solution, 100 mL
- Sodium thiosulfate solution, 100 mL

Nuclear fast red-aluminum sulfate solution 0.1 %

Cat. No.: 100121

Silver plating methods are highly specific. They are used to visualize fine structures such as microcalcification. Silver ions of the silver nitrate solution react with the carbonate and phosphate ions of the calcium and displace the calcium ions.




Components:

- Silver nitrate solution, 100 mL
- Sodium thiosulfate solution, 100 mL

Tb-color modified Kit for hot staining acid-fast bacteria

Cat. No.: 100497

Due to the high proportion of wax and lipids in the cell wall, mycobacteria absorb dyes "only very" slowly. To accelerate the uptake of fuchsin dye the carbolfuchsin solution applied to the specimen is normally heated until vapors evolve. The mycobacteria are thus referred to as being acid- and alcohol-fast in staining operations, and show up red on the microscope slide.




Components:

- Tb-color modified carbol-fuchsin solution, 500 mL
- Tb-color modified hydrochloric acid in ethanol, 2 x 500 mL
- Tb-color modified methylene blue solution, 500 mL

TB-fluor Kit for fluorescence staining acid-fast bacteria

Cat. No.: 109093

Fluorescence staining allows acid-fast rods to be easily identified even when their numbers are small. Depending on the filter, choose between two fluorescent dyes: Rhodamine B – red and Auramine O – yellowgreen. Potassium permanganate solution turns the background uniformly black in order to suppress troublesome intrinsic fluorescence of the specimen.



Components:

- Auramine-Rhodamine solution, 500 mL
- Decolorisation solution, 3 x 500 mL
- Counterstaining solution (KMnO₄), 2 x 500 mL

TISSUE-TROL™ Control Slides

Histology control slides are essential for formal quality assurance programs. Our slides are composed of paraffin-embedded tissue sections containing known characteristics and are used for monitoring staining performance in pathology laboratories.

Cat. No.	Product Description
TTR001	Acid Fast TISSUE-TROL™ Control Slides from mouse lung tissue containing <i>Mycobacterium gordonae</i>
TTR002	Amyloid TISSUE-TROL™ Control Slides from Human heart tissue containing intracellular amyloid
TTR013	Argentaffin TISSUE-TROL™ Control Slides from human intestine tissue containing argentaffin granules
TTR003	Elastic TISSUE-TROL™ Control Slides from human skin containing elastic fibers
TTR004	Fungi TISSUE-TROL™ Control Slides from mouse lung containing <i>Candida albicans</i>
TTR005	Gram Stain TISSUE-TROL™ Control Slides from mouse lung tissue containing <i>Staphylococcus aureus</i> and <i>Escherichia coli</i>
TTR006	Helicobacter TISSUE-TROL™ Control Slides from mouse intestine tissue containing <i>Helicobacter pylori</i>
TTR007	Iron TISSUE-TROL™ Control Slides from human liver containing intracellular and extracellular iron
TTR014	Melanin TISSUE-TROL™ Control Slides from human skin
TTR008	Mucin TISSUE-TROL™ Control Slides from human intestine tissue containing mucins
TTR009	PAS TISSUE-TROL™ Control Slides from human kidney
TTR010	Reticulum TISSUE-TROL™ Control Slides
TTR011	Spirochetes TISSUE-TROL™ Control Slides from rabbit testicle tissue containing spirochetes (<i>Treponema pallidum</i>)
TTR012	Trichrome TISSUE-TROL™ Control slides from human liver



Midas® III-Plus: The Cornerstone of Your Laboratory

The Midas® III-Plus Automated Stainer is a fully automated, benchtop stainer with a small footprint. It is designed to perform routine hematology and microbiology staining effortlessly. It utilizes the dip slide staining technique which ensures uniform staining for every sample. This reduces stain consumption up to 50%.

This versatile automated stainer combines reliable electronics and mechanics, now with a new optional plastic hood enclosure to make your laboratory a safer environment. Ideally suited for staining up to 300 slides per day, this stainer will become the new workhorse of your hematology/bacteriology lab.

Discover how you can get quality staining every time:

www.emdmillipore.com

Cat. No.	Product	Size	Package
64000H-94	Midas® III-Plus Slide Stainer	1EA	Chipboard Box
Additional Accessories			
64052E-94	Midas® Slide Carrier 30 slides	1EA	Chipboard Box
64054W-95	Midas® Vessels and Lids White 8/pk	1EA	Chipboard Box
64000WB-94	Midas® Water Bath	1EA	Chipboard Box
64000WBF-95	Midas® Water Bath Filters 2/pk	1EA	Chipboard Box
RW918003	Midas® 110V power cord (US)	1 EA	Chipboard Box

Immunohistochemistry



SNAP i.d.® 2.0 Protein Detection System

SNAP i.d.® 2.0 Protein Detection System for Immunohistochemistry (IHC) introduces a new capability to the innovative, vacuum-driven SNAP i.d.® 2.0 system. The IHC frame and slide holders allow you to block, probe and stain up to 12 tissue slides per frame. Reduced handling time and multiple-slide processing make this system ideal for antibody and protocol standardization.

Key Benefits

- Eliminates the need for pap pens
- Antibodies can be collected and reused
- Slide handling time is significantly decreased
- Less time spent on wash steps
- Parallel processing of multiple slides
- Efficient, vacuum-driven slide washing saves time and can reduce background

Base System

The SNAP i.d.® systems for IHC contain everything you need to get started, including the detection base, IHC frame and incubation cover, slide holders, an assembly fixture, vacuum tubing and a Quick Start User Guide.

Cat. No.	Product Description
SNAP2IHC	SNAP i.d.® 2.0 Protein Detection System - 1 IHC Frame (up to 12 slides)
SNAP2IHC2	SNAP i.d.® 2.0 Protein Detection System - 2 IHC Frames (up to 24 slides)

SNAP i.d.® 2.0 consumables

Cat. No.	Product Description
SNAP2FRIHC	SNAP i.d.® 2.0 IHC Frame
SNAP2SH	SNAP i.d.® 2.0 IHC Slide Holders 24/pk

Chromogenic Substrates

Chromogenic substrates and enhancers are used to localize and increase the staining intensity of tissue specimens stained using immunohistochemistry. Typically the enzymes horseradish peroxidase (HRP) and alkaline phosphatase (AP) are used in conjunction with the appropriate chromogen/substrate to produce a reaction product that is then visualized using light microscopy.

Cat. No.	Product Description
D0426	SIGMAFAST™ DAB with Metal Enhancer – SIGMAFAST™ DAB (3,3'-Diaminobenzidine tetrahydrochloride) with Metal Enhancer Tablets have been developed for use in immunohistology as a precipitating substrate for the localization of peroxidase activity. The DAB reaction is enhanced by the addition of CoCl ₂ and produces a more intense stain which is resistant to alcohol.
D3939	3,3'-Diaminobenzidine (DAB) Enhanced Liquid Substrate System tetrahydrochloride – The DAB (3,3'-Diaminobenzidine tetrahydrochloride) Two-component Enhanced Liquid Substrate System is a precipitating substrate system for peroxidase. DAB Two-component Liquid Substrate System produces a brown precipitate for use in immunohistological and immunoblotting staining procedures.
AEC101	AEC Staining Kit – The AEC (3-Amino-9-ethylcarbazole) Staining Kit contains universal reagents for staining peroxidase labeled compounds in immunohistochemistry and immunoblotting techniques. AEC produces an insoluble end product that is red in color.
F4648 F4523	SIGMAFAST™ Fast Red TR/Naphthol AS-MX Tablets – Developed for use as a precipitating substrate for the detection of alkaline phosphatase activity in immunohistology. Contains 0.6 mM levamisole to block endogenous alkaline phosphatase activity.
R3783	PermaRed/AP Kit – PermaRed/AP is a chromogenic substrate which produces a strong red color and is used for immunohistochemical (IHC) or <i>in situ</i> hybridization (ISH) detection of secondary antibodies conjugated to alkaline phosphatase (AP). After the enzymatic reaction is completed, the resulting dye is insoluble in organic solvents including alcohol and xylene, which allows tissue sections stained with PermaRed/AP to be dehydrated, cleared and mounted after staining.
SRE0079-50G SRE0079-250G	3,3'-Diaminobenzidine tetrahydrochloride hydrate C ₁₂ H ₁₄ N ₄ · 4HCl · xH ₂ O – Elite quality level DAB produces a brown precipitate when oxidized by peroxidase. Used in immunohistological and immunoblotting staining procedures.

Ancillary Reagents

Cat. No.	Qty.	Product Description
938B-03	50 mL predilute	Diamond: Antibody Diluent is a buffered solution for the dilution of both polyclonal and monoclonal antibodies as well as for the preparation of negative control reagents in IHC. It is designed to minimize the non-specific reaction that may be caused by the antibody and encourage specific antigen-antibody binding. Diamond: Antibody Diluent can also be used to stabilize diluted antibodies when stored at 2-8°C.
938B-05	200 mL predilute	
935B-06	200 mL concentrate	Tris Buffered Immunohistochemistry Wash Buffer + Tween® 20 is a 20X concentrated solution that is employed to rinse reagents off slides and to provide a medium for short-term storage of immunohistochemistry specimens in between application of reagents. When diluted, the ready to use solution is a 50mM Tris buffered solution with a pH range of 7.6±0.1.
935B-09	1000 mL concentrate	
954D-20	50 mL each predilute amplifier and polymer detector	The HiDef Detection™ Polymer System is a visualization system that is ready-to-use in immunohistochemical protocols. This two-step system uses an indirect method resulting in an antibody-enzyme complex that universally detects mouse and rabbit primary antibodies.
954D-30	100 mL each predilute amplifier and polymer detector	
957D-50	150 tests	DAB Substrate Kit is intended for use in immunohistochemistry (IHC) staining protocols. This reagent is designed to detect target antigens in formalin-fixed, paraffin-embedded tissue sections when used in conjunction with primary antibodies and appropriate secondary antibody conjugates in the IHC staining process.
957D-20	500 tests	
957D-60	1000 tests	
957D-30	2000 tests	When in the presence of HRP polymer, DAB produces a brown precipitate that can be readily visualized using light microscopy.
957D-40	5000 tests	

Antigen Retrieval

Due to the cross-linking chemistry used in fixation it is often necessary to pre-treat a fixed tissue sample prior to immunostaining to make target epitopes accessible. We offer enzymes and buffers to accommodate both enzymatic and HIER (Heat Induced Epitope Retrieval) methods of antigen retrieval.

Cat. No.	Product Description
SRE0063	New - Tris-EDTA Buffer pH 9.0, Ready-to-use, Antigen Retriever, IVD - For use in Diagnostic applications for the heat induced epitope (or antigen) retrieval (HIER) of formalin-fixed, paraffin-embedded (FFPE) tissues prior to immunohistochemical (IHC) staining.
C9999	Citrate Buffer, pH 6.0, 10x, Antigen Retriever - Used as a heat-induced antigen retriever on formalin-fixed paraffin-embedded (FFPE) tissue sections prior to application of antibodies.
E1161	EDTA Buffer, pH 8.5, 10x, Antigen Retriever - Used as a heat-induced antigen retriever on formalin-fixed paraffin-embedded (FFPE) tissue sections prior to application of antibodies.
T6455	Tris-HCl Buffer, pH 10, 10x, Antigen Retriever - Used as a heat-induced antigen retriever on formalin-fixed paraffin-embedded (FFPE) tissue sections prior to application of antibodies.
H3292	Histo/Zyme, pH 7.2, Ready-to-use, Antigen Retriever - Histo/Zyme is a stabilized proteolytic enzyme in buffered solution, used to unmask immunoreactive sites altered by fixation and/or the embedding process. Treatment with a proteolytic enzyme in immunohistochemistry (IHC) is required to digest excess linkages and expose the antigen. The use of Histo/Zyme results in superior staining over traditionally used enzyme pretreatments, such as pepsin and trypsin. Moreover, Histo/Zyme requires only a brief incubation at room temperature. Performing antigen recovery with Histo/Zyme will allow you to dilute the primary antibody 5 to 10 fold further and achieve optimal staining.
R2283	Pepsin - Used for proteolytic digestion of formalin-fixed paraffin-embedded (FFPE) tissue sections prior to application of antibodies.
920P	Trilogy™ IHC Pretreatment Solution. Combines 3 pretreatment steps: deparaffinization, rehydration, and unmasking

Featured Product

Trilogy™ – 3 Steps in 1

Trilogy™ is a product that combines the three pretreatment steps: deparaffinization, rehydration, and unmasking in immunohistochemistry stains. This product is completely biodegradable and non-toxic.

Trilogy™ may be used whenever immunohistochemistry staining is performed on formalin-fixed paraffin-embedded tissue sections. Using this product encourages standardization of the pretreatment procedure, thereby producing more consistent, more reliable results.



Research Antibodies


Analyzing protein expression in cells is imperative to understanding the cause and prevention of human diseases. Our comprehensive portfolio contains over 90,000 highly characterized primary antibodies plus nearly 2,500 secondary antibodies and conjugates. Over 37,000 of these antibodies have been validated for use in IHC applications. Our antibody offering is supported by scientific research—with more than 10,000 peer-reviewed citations per year—and our Antibody Guarantee, giving you confidence that we can support and drive your research from start to finish.

For Research Use Only. Not for use in diagnostic procedures.


Prestige Antibodies®

Our catalog of over 20,000 Prestige Antibodies® provides industry-leading depth of characterization. These exceptional antibodies were developed based on the goal of the Human Protein Atlas project to map the entire human tissue proteome via IHC staining. Each antibody is supported by over 500 IHC, IF, and/or Western blot images and all supporting data is publicly available on the HPA website. Thus, you can review countless images from a wide range of tissues to enable you to make confident purchasing decisions.


Cat. No.	Product Description
HPA005785	Anti-CD44 antibody produced in rabbit
HPA003259	Anti-MITF antibody produced in rabbit
HPA001328	Anti-BRAF antibody produced in rabbit
HPA003456	Anti-ZEB2 antibody produced in rabbit
HPA008405	Anti-COL1A1 antibody produced in rabbit
HPA042785	Anti-CSPG4 antibody produced in rabbit
HPA018530	Anti-EGFR antibody produced in rabbit
HPA001762	Anti-VIM antibody produced in rabbit
HPA001238	Anti-MMP9 antibody produced in rabbit
AMAB90627	Monoclonal Anti-HER2 antibody produced in mouse



Anti-CD44 (Product No. HPA005785): Immunohistochemical staining of human esophagus shows strong cytoplasmic and membranous positivity in squamous epithelial



Anti-BRAF (Product No. HPA001328) Immunohistochemical staining of human testis shows moderate cytoplasmic positivity in cells in seminiferous ducts.




Monoclonal Anti-HER2 (Product No. AMAB90627) Immunohistochemical staining of human HER2-positive breast cancer shows strong membranous positivity in tumor cells.

Discover more at


[SigmaAldrich.com/antibodies](https://www.sigmaaldrich.com/antibodies)

Cancer Research Antibodies

Recognizing both the tremendous opportunities and the challenges facing cancer researchers, we are dedicated to developing and refining tools and technologies for your study of cancer. We provide a comprehensive portfolio of cancer research antibodies extensively characterized and validated to support your application needs.



Anti-TWIST1 (Product No. SAB2108515): Immunohistochemistry of TWIST1 in Uterus tissue with TWIST1 antibody at 4-8 µg/mL.



Amyloid (β-A4-protein, PN: MAB348) staining on Alzheimer's Disease-Hypothalamus. Tissue pretreated with Citrate, pH 6.0. The antibody was diluted to 1:80, using IHC-Select® Detection with HRP-DAB. Immunoreactivity is seen as cerebral vascular amyloid deposition in the walls of blood vessels.

Cat. No.	Product Description
Signalling	
SAB4504331	Anti-phospho-Akt (pSer473) antibody
MAB5352	Anti-Notch 1 Antibody, clone mN1A
05-361	Anti-NFκB p52 Antibody
P4235	Anti-p53 antibody
05-516	Anti-Ras Antibody, clone RAS10
WH0003845M1	Monoclonal Anti-KRAS antibody
Growth Factors	
E2760	Monoclonal Anti-Epidermal Growth Factor Receptor antibody
E3138	Monoclonal Anti-Epidermal Growth Factor Receptor
V4262	Monoclonal Anti-Vascular Endothelial Growth Factor Receptor-1 antibody
V4758	Monoclonal Anti-Vascular Endothelial Growth Factor antibody
F6162	Anti-Fibroblast Growth Factor-Basic antibody
F5421	Anti-Fibroblast Growth Factor Receptor-1 antibody
Cell Death	
HPA045168	Anti-PARP1 antibody
A2856	Anti-ATG7 antibody
C8487	Anti-Caspase 3, Active antibody
C1354	Monoclonal Anti-Caspase 11 antibody
A8604	Monoclonal Anti-Annexin V antibody
Metabolism/Oxidative Stress	
T2949	Anti-mTOR antibody
U6382	Anti-UCP-1 Antibody
SAB4500003	Anti-BCL-2 antibody
G7781	Anti-Glutathione-S-Transferase (GST) antibody
S2147	Monoclonal Anti-Superoxide Dismutase (SOD) antibody

Neurodegenerative Pathology

The progressive loss of neuronal function can lead to serious debilitating diseases (e.g. Alzheimer's) and cause premature death. Count on our comprehensive portfolio of neurodegenerative antibodies so you can better explore and understand the mechanisms of neurological diseases.

Cat. No.	Product Description
Neuronal Cell Structure	
MAB377	Anti-NeuN Antibody, clone A60
T8578	Monoclonal Anti-β-Tubulin III (neuronal) antibody
MAB3420	Anti-Tau-1 Antibody, clone PC1C6
Neuronal Function	
AB152	Anti-Tyrosine Hydroxylase Antibody
AB144P	Anti-Choline Acetyltransferase Antibody
MAB5406	Anti-GAD67 Antibody, clone 1G10.2
Neuronal Development	
AB9610	Anti-Olig-2 Antibody
ABD69	Anti-Nestin, Human Antibody
AB1554	Anti-Nerve Growth Factor Receptor Antibody, p75
Neural Degeneration	
AB5320	Anti-NG2 Chondroitin Sulfate Proteoglycan Antibody
AB5078P	Anti-Beta-Amyloid 1-42 Antibody
AB5038	Anti-Synuclein α Antibody


Cell Marque™ IVD

Cell Marque™ IVD primary antibodies form a comprehensive portfolio of both novel and routine antibodies. Our team of experts includes medical professionals, which allows us to deliver the highest quality reagents and vast technical expertise to the IVD market.


Discover the entire product portfolio at:
[SigmaAldrich.com/tissuedx](https://www.sigmaaldrich.com/tissuedx)



Hematopathology



c-Myc (EP121)




SOX-11 (MRQ-58)


Cat. No.	Product Description
395R-18	c-Myc (EP121), 7 mL ready-to-use ready-to-use**
104R-18	CD4 (SP35), 7 mL ready-to-use**
138M-18	CD138/syndecan-1 (B-A38), 7 mL ready-to-use**
382M-18	SOX-11 (MRQ-58), 7 mL ready-to-use**
312R-18	PAX-5 (SP34), 7 mL ready-to-use**



Gastrointestinal (GI) Pathology



CDX-2 (EPR2764Y)




SATB2 (EP281)


Cat. No.	Product Description
235R-18	CDX-2 (EPR2764Y) 7 mL ready-to-use**
384R-18	SATB2 (EP281) 7 mL ready-to-use**



Dermatopathology



MART-1 (Melan A) (M2-7C10)




SOX-10 (EP268)


Cat. No.	Product Description
198M-18	CD123 (6H6) 7 mL ready-to-use**
281M-98	MART-1 (Melan A) (M2-7C10) 7 mL ready-to-use**
284M-98	Microphthalmia Transcription Factor (MiTF) (C5/D5) 7 mL ready-to-use**
383R-18	SOX-10 (EP268) 7 mL ready-to-use**



Genitourinary (GU) Pathology



ERG (EP111)




GATA3 (L50-823)


Cat. No.	Product Description
434R-18	ERG (EP111) 7 mL ready-to-use**
390M-18	GATA3 (L50-823) 7 mL ready-to-use**
441R-18	NKX3.1 (EP353) 7 mL ready-to-use**



Soft Tissue Pathology



MyoD1 (EP212)




TLE1 (1F5)


Cat. No.	Product Description
386R-18	MyoD1 (EP212) 7 mL ready-to-use**
354R-18	TFE3 (MRQ-37) 7 mL ready-to-use**
401M-18	TLE1 (1F5) 7 mL ready-to-use**



Breast/Gynecological Pathology



Ki-67 (SP6)



p120 Catenin (MRQ-5)

Cat. No.	Product Description
246R-18	E-cadherin (EP700Y) 7 mL ready-to-use**
275R-18	Ki-67 (SP6) 7 mL ready-to-use**
420M-18	p120 Catenin (MRQ-5) 7 mL ready-to-use**

**Also available in concentrate sizes.



Water Purification Solutions

Our Milli-Q® lab water portfolio offers a broad range of pure and ultrapure water purification systems and services designed for scientists working on tissue pathology. Our pure, ultrapure, and Clinical Laboratory Reagent Water (CLRW) purification systems deliver constant and reliable water quality in the volumes and distribution formats needed to best fit your pathology laboratory's specific needs.

- Milli-Q® IQ 7000 purification system provides ultrapure water free of trace organic molecules, including enzymes and microbials, and ions.
- Elix® and Milli-Q® HX 7000 purification systems contain Elix® electrodeionization technology which ensures consistent and reliable pure water, free of contaminants that may interfere with your fixation and staining procedures.
- Milli-Q® CLX 7000 series analyzer feed systems containing Elix® technology can feed automated, high-volume, high-throughput slide stainers, and can provide CLRW to your laboratory's clinical chemistry analyzer.

To learn more about our portfolio of biomedical lab water systems, please visit:
[SigmaAldrich.com/clinical-water](https://www.sigmaaldrich.com/clinical-water)

Products included in this brochure are not available in all countries.
Contact your local sales representative or dealer for details.

MilliporeSigma
400 Summit Drive
Burlington, MA 01803

To place an order or receive technical assistance

Order/Customer Service: SigmaAldrich.com/order
Technical Service: SigmaAldrich.com/techservice
Safety-related Information: SigmaAldrich.com/safetycenter

SigmaAldrich.com

