





MILLIPORE  
SIGMA

# BENCH TO BULK

Polyethylene Glycol (PEG)  
and Polyethylene Oxide (PEO)  
for Medical Devices



MilliporeSigma is the U.S.  
and Canada Life Science  
business of Merck KGaA,  
Darmstadt, Germany.

**Sigma-Aldrich**<sup>®</sup>  
Lab & Production Materials

# Polyethylene Glycol (PEG) and Polyethylene Oxide (PEO) for Medical Devices

Due to their solubility, biocompatibility, non-toxicity, and low immunogenicity, PEG and PEO materials are used in a wide variety of medical devices and medical device applications (Page 4)

## Let us support you with:

- Comprehensive portfolio of over 450 well-defined PEGs available in quantities that support a seamless scale-up from R&D to production
- Wide range of molecular weights, end functionalities, reactivities, and polymer architecture for medical device development and production
- Consistent, reliable and transparent supply chain for risk mitigation, including dual-source and inventory management. For details, visit [SigmaAldrich.com/mclarity](https://www.sigmaaldrich.com/mclarity)
- Custom packaging, additional QC testing and potential customization of PEG chemistry

## Available product types:

- Heterobifunctional PEGs
- Homobifunctional PEGs
- Monofunctional PEGs
- PEG Dendrimers and Multi-Arm PEGS
- PEG Copolymers
- PEG Analytical Standards for GPC
- PEG and Oligo Ethylene Glycol
- Poly(ethylene oxide) (PEO) Powders
- High Oligomer Purity PEGs

### Functionality

- **Monofunctional:** PEGylation, surface conjugation and nanoparticle coating
- **Homo/heterofunctional:** Conjugation & crosslinking for hydrogels

### Polymer Architecture

- **Linear:** PEGylation, bioconjugation and crosslinking
- **Multiarm (4-, 6-, 8-arm):** Crosslinked into hydrogels and scaffolds for drug delivery or tissue engineering
- **Y-shaped:** PEGylation, as the branched structure may improve stability *in vivo*

### Reactivity

- **Reactive end groups, such as an N-hydroxysuccinimide ester, thiol, or carboxyl groups:** Covalent conjugation
- **Azide or alkyne reactive groups:** Click chemistry for conjugation or hydrogen formation
- **Acrylate-terminate:** Polymerization and photopolymerization under mild reactive conditions

### Molecular Weight

- **$M_w \leq \text{kDa}$ :** PEGylation of proteins
- **$M_w \geq \text{kDa}$ :** Hydrogel formation of typical mesh size and mechanical properties, conjugation to small molecules, siRNA, and peptides
- **$M_w < \text{kDa}$ :** Surface conjugation and crosslinking

## HETEROBIFUNCTIONAL PEGS

Cat. No.	Name	Average Mol. Wt./Purity
407348	Poly(ethylene glycol) phenyl ether acrylate	average M <sub>n</sub> 324, contains 150-350 ppm MEHQ as inhibitor
409537	Poly(ethylene glycol) methacrylate	average M <sub>n</sub> 360, contains 500-800 ppm MEHQ as inhibitor
757837	Poly(ethylene glycol) 2-mercaptoethyl ether acetic acid	average M <sub>n</sub> 3,500
757810	Poly(ethylene glycol) 2-mercaptoethyl ether acetic acid	average M <sub>n</sub> 1,000
757896	Poly(ethylene glycol) 2-aminoethyl ether acetic acid	average M <sub>n</sub> 3,500
757799	Poly(ethylene glycol) (N-hydroxysuccinimide 5-pentanoate) ether 2-(biotinylamino)ethane	average M <sub>n</sub> 3,800

Cat. No.	Name	Average Mol. Wt./Purity
757756	Poly(ethylene glycol) 2-aminoethyl ether biotin	average M <sub>n</sub> 2,300
409529	Poly(ethylene glycol) methacrylate	average M <sub>n</sub> 500, contains 900 ppm monomethyl ether hydroquinone as inhibitor
757845	Poly(ethylene glycol) 2-mercaptoethyl ether acetic acid	PEG average M <sub>n</sub> 5,000 (n~110)
757888	Poly(ethylene glycol) 2-aminoethyl ether acetic acid	average M <sub>n</sub> 2,100
757829	Poly(ethylene glycol) 2-mercaptoethyl ether acetic acid	average M <sub>n</sub> 2,100
757918	Poly(ethylene glycol) 2-aminoethyl ether acetic acid	average M <sub>n</sub> 5,000

## HOMOBIFUNCTIONAL PEGS

Cat. No.	Name	Average Mol. Wt./Purity
407038	Poly(ethylene glycol) bis(carboxymethyl) ether	average M <sub>n</sub> 600
406996	Poly(ethylene glycol) bis(carboxymethyl) ether	average M <sub>n</sub> 250
335681	Ethylene glycol dimethacrylate	98%, contains 90-110 ppm monomethyl ether hydroquinone as inhibitor
409510	Poly(ethylene glycol) dimethacrylate	average M <sub>n</sub> 550, contains 80-120 ppm MEHQ as inhibitor, 270-330 ppm BHT as inhibitor
475696	Poly(ethylene glycol) diglycidyl ether	average M <sub>n</sub> 500
445886	Poly(ethylene glycol) dimethyl ether	average M <sub>n</sub> ~500, contains 100 ppm BHT as stabilizer
398802	Tetra(ethylene glycol) diacrylate	technical grade, contains 150-200 ppm MEHQ as inhibitor, 100-150 ppm HQ as inhibitor
305413	Poly(ethylene glycol) distearate	average M <sub>n</sub> ~930
701963	Poly(ethylene glycol) diacrylate	average M <sub>n</sub> 6,000, contains ≤1500 ppm MEHQ as inhibitor
14526	O,O'-Bis(2-aminopropyl) polypropylene glycol-block-polyethylene glycol-block-polypropylene glycol	500

Cat. No.	Name	Average Mol. Wt./Purity
458074	Poly(ethylene glycol) dibenzoate	average M <sub>n</sub> ~410
753084	Poly(ethylene glycol) diamine	average M <sub>n</sub> 2,000
452572	Poly(ethylene glycol) bis (3-aminopropyl) terminated	M <sub>n</sub> ~1,500
261548	Triethylene glycol dimethacrylate	contains 80-120 ppm MEHQ as inhibitor, 95%
729094	Poly(ethylene glycol) diacrylate	average M <sub>n</sub> 10,000, contains MEHQ as inhibitor
329800	Tri(ethylene glycol) divinyl ether	98%
805505	Poly(ethylene glycol) diglycidyl ether	M <sub>n</sub> 1,000
437433	Di(ethylene glycol) diacrylate	technical grade, 75%
14529	O,O'-Bis(2-aminopropyl) polypropylene glycol-block-polyethylene glycol-block-polypropylene glycol	1,900
437468	Poly(ethylene glycol) dimethacrylate	average M <sub>n</sub> 750, contains 900-1100 ppm MEHQ as inhibitor
172405	Tetraethylene glycol dimethyl ether	≥ 99%

## MONOFUNCTIONAL PEGS

Cat. No.	Name	Average Mol. Wt./Purity
457876	Poly(ethylene glycol) methyl ether methacrylate solution	average M <sub>n</sub> 2,000, 50 wt. % in H <sub>2</sub> O
202509	Poly(ethylene glycol) methyl ether	average M <sub>n</sub> ~2,000
202495	Poly(ethylene glycol) methyl ether	average M <sub>n</sub> 750
90450	Triethylene glycol monomethyl ether	purum, ≥97.0% (GC)
81323	Poly(ethylene glycol) methyl ether	average M <sub>n</sub> 5,000
729159	Poly(ethylene glycol) methyl ether thiol	average M <sub>n</sub> 6,000
732613	Poly(ethylene glycol) methyl ether	average M <sub>n</sub> 20,000

Cat. No.	Name	Average Mol. Wt./Purity
729140	Poly(ethylene glycol) methyl ether thiol	average M <sub>n</sub> 2,000
202487	Poly(ethylene glycol) methyl ether	average M <sub>n</sub> 550
317292	Triethylene glycol monomethyl ether	95%
732621	Poly(ethylene glycol) methyl ether	average M <sub>n</sub> 10,000
447943	Poly(ethylene glycol) methyl ether methacrylate	average M <sub>n</sub> 500, contains 200 ppm BHT as inhibitor, 100 ppm MEHQ as inhibitor
729108	Poly(ethylene glycol) methyl ether thiol	average M <sub>n</sub> 800
447935	Poly(ethylene glycol) methyl ether methacrylate	average M <sub>n</sub> 300, contains 100 ppm MEHQ as inhibitor, 300 ppm BHT as inhibitor

## MULTIARM PEGS

Cat. No.	Name	Average Mol. Wt./Purity
760919	Poly(ethylene glycol), 8 acetylene dendron, generation 2	average M <sub>n</sub> 6,900

Cat. No.	Name	Average Mol. Wt./Purity
441864	Glycerol ethoxylate	average M <sub>n</sub> ~1,000

## PEG & OLIGO(ETHYLENE GLYCOL)

Cat. No.	Name	Average Mol. Wt./Purity
81260	Poly(ethylene glycol)	average M <sub>n</sub> 6,000
182028	Poly(ethylene oxide)	average M <sub>v</sub> 600,000, powder
P2139	Poly(ethylene glycol)	average mol M <sub>w</sub> 8,000, powder
110175	Tetraethylene glycol	99%
181986	Poly(ethylene oxide)	average M <sub>v</sub> 100,000, powder
P4338	Poly(ethylene glycol)	BioXtra, average M <sub>w</sub> 3,350, powder
189456	Poly(ethylene oxide)	average M <sub>v</sub> ~900,000, powder
202401	Poly(ethylene glycol)	average M <sub>n</sub> 600, waxy solid (moist)
372781	Poly(ethylene oxide)	average M <sub>v</sub> ~1,000,000, powder

Cat. No.	Name	Average Mol. Wt./Purity
295906	Poly(ethylene glycol)	average M <sub>n</sub> 2,050, chips
81310	Poly(ethylene glycol)	M <sub>w</sub> 35,000
81240	Poly(ethylene glycol)	average M <sub>n</sub> 4,000, platelets
202371	Poly(ethylene glycol)	average M <sub>n</sub> 300
181994	Poly(ethylene oxide)	average M <sub>v</sub> 200,000 (nominal), powder
202398	Poly(ethylene glycol)	average M <sub>n</sub> 400
81300	Poly(ethylene glycol)	average M <sub>n</sub> 20,000
202444	Poly(ethylene glycol)	average M <sub>n</sub> 3,350, powder
81210	Poly(ethylene glycol)	average M <sub>w</sub> 1,500
P3265	Poly(ethylene glycol)	average M <sub>w</sub> 400
202436	Poly(ethylene glycol)	average M <sub>n</sub> 1,305-1,595, waxy solid

## Key Industrial Applications for PEGs

### Medical Devices, Dental Implants, Biosensors, Diagnostics Applications

Cat. No.	Name	Average Mol. Wt./Purity
189456	Poly(ethylene oxide)	average M <sub>v</sub> ~900,000, powder
81310	Poly(ethylene glycol)	M <sub>w</sub> 35,000
181986	Poly(ethylene oxide)	average M <sub>v</sub> 100,000, powder
81210	Poly(ethylene glycol)	average M <sub>w</sub> 1,500
P2139	Poly(ethylene glycol)	average mol M <sub>w</sub> 8,000, powder
261548	Triethylene glycol dimethacrylate	contains 80-120 ppm MEHQ as inhibitor, 95%

Cat. No.	Name	Average Mol. Wt./Purity
202401	Poly(ethylene glycol)	average M <sub>n</sub> 600, waxy solid (moist)
372781	Poly(ethylene oxide)	average M <sub>v</sub> ~1,000,000, powder
445886	Poly(ethylene glycol) dimethyl ether	average M <sub>n</sub> ~500, contains 100 ppm BHT as stabilizer
202371	Poly(ethylene glycol)	average M <sub>n</sub> 300
P4338	Poly(ethylene glycol)	BioXtra, average mol M <sub>w</sub> 3,350, powder

### Ophthalmic Contact Lens Applications

Cat. No.	Name	Average Mol. Wt./Purity
407348	Poly(ethylene glycol) phenyl ether acrylate	average M <sub>n</sub> 324, contains 150-350 ppm MEHQ as inhibitor
335681	Ethylene glycol dimethacrylate	98%, contains 90-110 ppm monomethyl ether hydroquinone as inhibitor

Cat. No.	Name	Average Mol. Wt./Purity
398802	Tetra(ethylene glycol) diacrylate	technical grade, contains 150-200 ppm MEHQ as inhibitor, 100-150 ppm HQ as inhibitor
437433	Di(ethylene glycol) diacrylate	technical grade, 75%

### Filtration Membrane Applications

Cat. No.	Name	Average Mol. Wt./Purity
409537	Poly(ethylene glycol) methacrylate	average M <sub>n</sub> 360, contains 500-800 ppm MEHQ as inhibitor

Cat. No.	Name	Average Mol. Wt./Purity
81260	Poly(ethylene glycol)	average M <sub>n</sub> 6,000
475696	Poly(ethylene glycol) diglycidyl ether	average M <sub>n</sub> 500

### Optical Coating Applications

Cat. No.	Name	Average Mol. Wt./Purity
202371	Poly(ethylene glycol)	average M <sub>n</sub> 300

Cat. No.	Name	Average Mol. Wt./Purity
90450	Triethylene glycol monomethyl ether	purum, ≥97.0% (GC)
202398	Poly(ethylene glycol)	average M <sub>n</sub> 400

### Hydrogel Applications

Cat. No.	Name	Average Mol. Wt./Purity
729086	Poly(ethylene glycol) diacrylate	average M <sub>n</sub> 1,000, contains MEHQ as inhibitor
701963	Poly(ethylene glycol) diacrylate	average M <sub>n</sub> 6,000, contains ≤1500 ppm MEHQ as inhibitor
701963	Poly(ethylene glycol) diacrylate	average M <sub>n</sub> 6,000, contains ≤1500 ppm MEHQ as inhibitor
729094	Poly(ethylene glycol) diacrylate	average M <sub>n</sub> 10,000, contains MEHQ as inhibitor

Cat. No.	Name	Average Mol. Wt./Purity
767549	Poly(ethylene glycol) diacrylate	average M <sub>n</sub> 20,000, contains ≤1000 ppm MEHQ as inhibitor
409510	Poly(ethylene glycol) dimethacrylate	average M <sub>n</sub> 550, contains 80-120 ppm MEHQ as inhibitor, 270-330 ppm BHT as inhibitor
202495	Poly(ethylene glycol) methyl ether	average M <sub>n</sub> 750

# Sigma-Aldrich®

Lab & Production Materials

MilliporeSigma  
400 Summit Drive  
Burlington, MA 01803

Visit [SigmaAldrich.com/peg](https://SigmaAldrich.com/peg) for more PEG products and resources

Visit [SigmaAldrich.com/pegselection](https://SigmaAldrich.com/pegselection) to read technical content to aid in your PEG selection and knowledge

## To place an order or receive technical assistance

Order/Customer Service: [SigmaAldrich.com/order](https://SigmaAldrich.com/order)

Technical Service: [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice)

Safety-related Information: [SigmaAldrich.com/safetycenter](https://SigmaAldrich.com/safetycenter)

[SigmaAldrich.com](https://SigmaAldrich.com)

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. MilliporeSigma, the vibrant M, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MS\_BR9359EN Ver. 1.0

41179

05/2022